



碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-168-3301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

Alexa Fluor 488标记山羊抗小鼠IgG (H+L)

| 产品编号 | 产品名称 | 包装 |
|-------|---------------------------------|-------|
| A0428 | Alexa Fluor 488标记山羊抗小鼠IgG (H+L) | 100μl |

产品简介:

- 本Alexa Fluor 488标记山羊抗小鼠IgG (H+L) (Alexa Fluor 488-labeled Goat Anti-Mouse IgG (H+L))为进口分装, 用于免疫荧光染色。
- Alexa Fluor 488是一种常用的非常明亮的绿色荧光探针。它比绝大部分常用的绿色荧光探针更加明亮, 更加不容易淬灭, 而且背景更低。Alexa Fluor 488的荧光光谱与FITC和Cy2比较接近。Alexa Fluor 488的吸收(激发)和发射峰参见下表。

| Fluorophore | Absorption Peak (nm) | Emission Peak (nm) |
|-----------------|----------------------|--------------------|
| Alexa Fluor 488 | 495 | 519 |

- 本抗体为用纯化的小鼠IgG免疫山羊, 然后用亲和纯化柱对获得的抗血清进行纯化, 并经过人IgG、马IgG、牛IgG(bovine IgG)、兔IgG和猪IgG吸附纯化的优质二抗。对人IgG、马IgG、牛IgG(bovine IgG)、兔IgG和猪IgG几乎没有结合能力。特别适合于对于二抗种属特异性要求比较高的荧光染色实验。
- 本Alexa Fluor 488标记山羊抗小鼠IgG (H+L)用于免疫荧光染色时的推荐稀释比例为1:500。实际实验操作过程中需根据抗原和抗体的具体情况适当调节荧光标记二抗的稀释比例, 推荐的调节范围为1:200-1000。
- 本抗体如果用于常规的免疫染色, 以每次检测需1毫升1:500稀释的荧光标记二抗计, 至少可以检测50次。如果适当重复使用已经使用过的荧光标记二抗, 至少可以多检测150-250次。

包装清单:

| 产品编号 | 产品名称 | 包装 |
|-------|---------------------------------|-------|
| A0428 | Alexa Fluor 488标记山羊抗小鼠IgG (H+L) | 100μl |
| — | 说明书 | 1份 |

保存条件:

-20°C避光保存, 一年有效。

注意事项:

- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明:

1. 免疫荧光染色请参考相关实验步骤进行。起始稀释浓度按照产品简介中推荐的稀释比例进行稀释。
2. 如果希望重复使用稀释的荧光标记二抗, 稀释的荧光标记二抗4°C保存。

使用本产品的文献:

1. 刘军楼, 汪悦, 徐力, 杨继兵, 孙志岭, 邓海山, 徐建亚. 冬凌草甲素对人胰腺癌细胞增殖的抑制作用及其机制. 世界华人消化杂志. 2010年 11月;18(32):3407-3412.
2. Yang HB, Yang X, Cao J, Li S, Liu YN, Suo ZW, Cui HB, Guo Z, Hu XD. cAMP-dependent protein kinase activated Fyn in spinal dorsal horn to regulate NMDA receptor function during inflammatory pain. J Neurochem. 2011 Jan;116(1):93-104.
3. Peng H, Wu Y, Zhang Y. Efficient delivery of DNA and morpholinos into mouse preimplantation embryos by electroporation. PLoS One. 2012; 7(8):e43748.
4. Wang X, Bai J, Zhang L, Wang X, Li Y, Jiang P. Poly(A)-binding protein interacts with the nucleocapsid protein of porcine reproductive and respiratory syndrome virus and participates in viral replication. Antiviral Res. 2012 Dec;96(3):315-23.
5. Li X, Li Z, Li N, Qi J, Fan K, Yin P, Zhao C, Liu Y, Yao W, Cai X, Wang L, Zha X. MAGI2 enhances the sensitivity of BEL-7404 human hepatocellular carcinoma cells to staurosporine-induced apoptosis by increasing PTEN stability. Int J Mol Med. 2013 Aug;32(2):439-47.
6. Wang HQ, Meng S, Li ZR, Peng ZG, Han YX, Guo SS, Cui XL, Li YH, Jiang JD. The antiviral effect of 7-hydroxyisoflavone against Enterovirus 71 in vitro. J Asian Nat Prod Res. 2013 Apr;15(4):382-9.
7. Tu C, Li N, Zhu L, Zhou L, Su Y, Li P, Zhu X. Cationic long-chain hyperbranched poly(ethylene glycol)s with low charge density for gene delivery. Polym. Chem. 2013;4:393-401.
8. Yang Z, Liu J, Liu H, Qiu M, Liu Q, Zheng L, Pang M, Quan F, Zhang Y. Isolation and Characterization of SSEA3+ Stem Cells Derived from Goat Skin Fibroblasts. Cell Reprogram. 2013 Jun;15(3):195-205.
9. Cai QQ, Zheng LF, Fan RF, Lian H, Zhou L, Song HY, Tang YY, Feng XY, Guo ZK, Wang ZY, Zhu JX. Distribution of dopamine receptors D1- and D2-immunoreactive neurons in the dorsal motor nucleus of vagus in rats. Auton Neurosci. 2013 Jun;176(1-2):48-53.
10. Li D, Liu Y, Yang Y, Chen JH, Yang J, Zou LY, Tian ZQ, Lv J, Xia PY. Looped Host Defense Peptide CLP-19 Binds to Microtubules and Inhibits Surface Expression of TLR4 on Mouse Macrophages. J Immunol. 2013 Jun 15;190(12):6083-92.
11. Gao C, Ding Y, Zhong L, Jiang L, Geng C, Yao X, Cao J. Tacrine induces apoptosis through lysosome- and mitochondria-dependent pathway in HepG2 cells. Toxicol In Vitro. 2014 Jun;28(4):667-74.

12. Jin Y, Sun C, Feng L, Li P, Xiao L, Ren Y, Wang D, Li C, Chen L. Regulation of SIV antigen-specific CD4⁺ T cellular immunity via autophagosome-mediated MHC II molecule-targeting antigen presentation in mice. *PLoS One*. 2014 Mar 26;9(3):e93143.
13. Mao J, Liu S. Regulation of RPE barrier function by all-trans retinoic acid in myopia. *Neurosci Lett*. 2014 May 7;568:17-22.
14. Miao SH, Sun HB, Ye Y, Yang JJ, Shi YW, Lu M, Hu G, Zhou JW. Astrocytic JWA Expression is Essential to Dopaminergic Neuron Survival in the Pathogenesis of Parkinson's Disease. *CNS Neurosci Ther*. 2014 Aug;20(8):754-62.
15. Wen J, Pei H, Wang X, Xie C, Li S, Huang L, Qiu N, Wang W, Cheng X, Chen L. 3-Gambogic acid exhibits anti-psoriatic efficacy through inhibition of angiogenesis and inflammation. *J Dermatol Sci*. 2014 Jun;74(3):242-50.
16. Liu FL, Kuang YQ, Mu D, Zheng HY, Zhu JW, Zheng YT. The effect of exon 7 deletion during the evolution of TRIMCyp fusion proteins on viral restriction, cytoplasmic body formation and multimerization. *PLoS One*. 2015 Mar 30;10(3):e0121666.
17. Chen H, Zhang L, Guo Z, Wang Y, He R, Qin Y, Quan F, Zhang Y. Improving the development of early bovine somatic-cell nuclear transfer embryos by treating adult donor cells with vitamin C. *Mol Reprod Dev*. 2015 Nov;82(11):867-79.
18. Yang X, Ndawula C Jr, Zhou H, Gong X, Jin J. JF-305, a pancreatic cancer cell line is highly sensitive to the PARP inhibitor olaparib. *Oncol Lett*. 2015 Feb;9(2):757-761.
19. Li H, Chen L, Zhang M, Zhang B. Foxa1 gene and protein in developing rat eccrine sweat glands. *J Mol Histol*. 2016 Oct 27. [Epub ahead of print]
20. Chen H, Zhang L, Deng T, Zou P, Wang Y, Quan F, Zhang Y. Effects of oocyte vitrification on epigenetic status in early bovine embryos. *Theriogenology*. 2016 Aug;86(3):868-78.
21. Liu YH, Jin JL, Wang YZ, Tan Y, Zhou YY, Peng T, Li F, Liang WD, Chartrand P, Jiang YY, Shen ZF. Protrusion-localized STAT3 mRNA promotes metastasis of highly metastatic hepatocellular carcinoma cells in vitro. *Acta Pharmacol Sin*. 2016 Jun;37(6):805-13.
22. Wang Y, Tu L, Li Y, Chen D, Wang S. Notoginsenoside R1 Protects against Neonatal Cerebral Hypoxic-Ischemic Injury through Estrogen Receptor-Dependent Activation of Endoplasmic Reticulum Stress Pathways. *J Pharmacol Exp Ther*. 2016 Jun;357(3):591-605.
23. Wang B, Wang D, Zhu Z, Wang W, Zhang X, Tang F, Zhou Y, Wang H, Liu M, Yao X, Yan X. The role of extracellular-5'-nucleotidase/CD73 in glioma peritumoural brain edema. *Neurol Sci*. 2016 Apr;37(4):603-11.
24. Chen C, Wang X, Huang X, Yong H, Shen J, Tang Q, Zhu J, Ni J, Feng Z. Nicotinamide N-methyltransferase: a potential biomarker for worse prognosis in gastric carcinoma. *Am J Cancer Res*. 2016 Feb 15;6(3):649-63.

Version 2016.10.09