



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-168-3301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

Cy3标记驴抗山羊IgG (H+L)

产品编号	产品名称	包装
A0502	Cy3标记驴抗山羊IgG (H+L)	50μl

产品简介:

- 本Cy3标记驴抗山羊IgG (H+L) (Cy3-labeled Donkey Anti-Goat IgG (H+L))为进口分装, 用于免疫荧光染色。
- Cy3是近年新开发的一种红色荧光探针。它比绝大部分其它的红色荧光探针更加明亮, 更加不容易淬灭, 而且背景更低。Cy3的吸收(激发)和发射峰参见下表。

Fluorophore	Absorption Peak (nm)	Emission Peak (nm)
Cy3	550	570

- 本抗体为用纯化的山羊IgG免疫驴, 然后用亲和纯化柱对获得的抗血清进行纯化, 并经过小鼠、大鼠、兔、牛和人等多种IgG反复吸附纯化的高品质二抗。对于经过吸附的这些物种的IgG几乎没有结合能力。
- 本Cy3标记驴抗山羊IgG (H+L)用于免疫荧光染色时的推荐稀释比例为1:500。实际实验操作过程中需根据抗原和抗体的具体情况适当调节荧光标记二抗的稀释比例, 推荐的调节范围为1:200-1000。
- 本抗体如果用于常规的免疫染色, 以每次检测需1毫升1:500稀释的荧光标记二抗计, 至少可以检测25次。如果适当重复使用已经使用过的荧光标记二抗, 至少可以多检测75-125次。

包装清单:

产品编号	产品名称	包装
A0502	Cy3标记驴抗山羊IgG (H+L)	50μl
—	说明书	1份

保存条件:

-20°C避光保存, 一年有效。

注意事项:

- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明:

1. 免疫荧光染色请参考相关实验步骤进行。起始稀释浓度按照产品简介中推荐的稀释比例进行稀释。
2. 如果希望重复使用稀释的荧光标记二抗, 稀释的荧光标记二抗4°C保存。

使用本产品的文献:

1. Cai H, Yin D, Zhang L, Yang X, Xu X, Liu W, Zheng X, Zhang H, Wang J, Xu Y, Cheng D, Zheng M, Han Y, Wu M, Wang Y. Preparation and biological evaluation of 2-amino-6-[18F]fluoro-9-(4-hydroxy-3-hydroxymethylbutyl) purine (6-[18F]FPCV) as a novel PET probe for imaging HSV1-tk reporter gene expression. Nucl Med Biol. 2007 Aug; 34(6):717-25.
2. Ye JX, Zhang YT, Zhang XG, Ren DM, Chen WC. Recombinant attenuated Salmonella harboring 4-1BB ligand gene enhances cellular immunity. Vaccine. 2009 Mar 10;27(11):1717-23.
3. Yang F, Yang YP, Mao CJ, Cao BY, Cai ZL, Shi JJ, Huang JZ, Zhang P, Liu CF. Role of autophagy and proteasome degradation pathways in apoptosis of PC12 cells overexpressing human alpha-synuclein. Neurosci Lett. 2009 May 1;454(3):203-8.
4. Huang JZ, Chen YZ, Su M, Zheng HF, Yang YP, Chen J, Liu CF. dl-3-n-Butylphthalide prevents oxidative damage and reduces mitochondrial dysfunction in an MPP(+)-induced cellular model of Parkinson's disease. Neurosci Lett. 2010;475(2):89-94.
5. 吴毅梅, 江山平, 吕志强, 张蔚, 戴冽, 郑东辉. 瞬时感受器电位M7通道对RBL-2H3细胞增殖及凋亡的影响. 中山大学学报: 医学科学版. 2010 Jun;40(3):338-342.
6. Gan SW, Ran JH, Chen H, Ren ZQ, Sun SQ, Zhu SJ, Lu WT, Xu J, Zhang B, Huang J, Wang KJ, Chen Z. Lysosomal degradation of retinal glial AQP4 following its internalization induced by acute ocularhypertension. Neurosci Lett. 2012 May 10;516(1):135-40.
7. Ma F, Li Y, Jia L, Han Y, Cheng J, Li H, Qi Y, Du J. Macrophage-stimulated cardiac fibroblast production of IL-6 is essential for TGF β/Smad activationand cardiac fibrosis induced by angiotensin II. PLoS One. 2012;7(5):e35144.
8. Wang Y, Wang K, Chao R, Li J, Zhou L, Ma J, Yan J. Neuroprotective effect of vaccination with autoantigen-pulsed dendritic cells after spinal cordinjury. J Surg Res. 2012 Jul;176(1):281-92.
9. Chen Y, Li ZY, Yang Y, Zhang HJ. Uncoupling protein 2 regulates glucagon-like peptide-1 secretion in L-cells. World J Gastroenterol. 2012 Jul 14;18(26):3451-7.
10. Liu L, Zhang C, Hu Y, Peng B. Protective effect of metformin on periapical lesions in rats by decreasing the ratio of receptor activator of nuclearfactor kappa B ligand/osteoprotegerin. J Endod. 2012 Jul;38(7):943-7.
11. Song N, Zhang G, Geng W, Liu Z, Jin W, Li L, Cao Y, Zhu D, Yu J, Shen L. Acid sensing ion channel 1 in lateral hypothalamus contributes to breathing control. PLoS One. 2012;7(7):e39982.
12. Lv J, Zhao Z, Chen Y, Wang Q, Tao Y, Yang L, Fan TP, Liu C. The chinese herbal decoction danggui buxue tang inhibits angiogenesis in a rat model of liver fibrosis. Evid Based Complement Alternat Med.

2012;2012:284963.

13. Liu AJ, Zang P, Guo JM, Wang W, Dong WZ, Guo W, Xiong ZG, Wang WZ, Su DF. Involvement of acetylcholine- α 7nAChR in the protective effects of arterial baroreflex against ischemic stroke. *CNS Neurosci Ther.* 2012 Nov;18(11):918-26.
14. Zhang SL, Zheng DJ, Fan WZ, Wei DX, Peng SW, Tang MM, Chen GQ, Wei CJ. Transient embolization with microspheres of polyhydroxyalkanoate renders efficient adenoviral transduction of pancreatic capillary in vivo. *J Gene Med.* 2012 Aug;14(8):530-9.
15. Liu W, Yin L, Dai Y. Modification of SR-PSOX functions by multi-point mutations of basic amino acid residues. *Biochimie.* 2013 Feb;95(2):224-30.
16. Wang L, Sun Z, Liu L, Peng B. Expression of CX3CL1 and its receptor, CX3CR1, in the development of periapical lesions. *Int Endod J.* 2013 May 29.
17. Yang X, Zhang Z, Wang X, Wang Y, Zhang X, Lu H, Wang SL. Cytochrome P450 2A13 enhances the sensitivity of human bronchial epithelial cells to aflatoxin B1-induced DNA damage. *Toxicol Appl Pharmacol.* 2013 Jul 15;270(2):114-21.
18. Luo D, Liu D, Zhou X, Yang S, Tang C, Liu G. Bone marrow-derived cells contribute to NDEA-induced lung squamous cell carcinoma. *Tumour Biol.* 2013 Feb;34(1):145-54.
19. Wan J, Zhou J, Zhao H, Wang M, Wei Z, Gao H, Wang Y, Cui H. Sonic hedgehog pathway contributes to gastric cancer cell growth and proliferation. *Biores Open Access.* 2014 Apr 1;3(2):53-9.
20. Li Y, Wang YQ, Wang DH, Hou WP, Zhang Y, Li M, Li FR, Mu J, Du X, Pang F, Yuan FH. Costimulatory molecule VSIG4 exclusively expressed on macrophages alleviates renal tubulointerstitial injury in VSIG4 KO mice. *J Nephrol.* 2014 Feb;27(1):29-36.
21. Li H, Chen L, Zeng S, Li X, Zhang X, Lin C, Zhang M, Xie S, He Y, Shu S, Yang L, Tang S, Fu X. Matrigel basement membrane matrix induces eccrine sweat gland cells to reconstitute sweat gland-like structures in nude mice. *Exp Cell Res.* 2015 Mar 1;332(1):67-77.
22. Wang JJ, Liu YL, Sun YC, Ge W, Wang YY, Dyce PW, Hou R, Shen W. Basic Fibroblast Growth Factor Stimulates the Proliferation of Bone Marrow Mesenchymal Stem Cells in Giant Panda (*Ailuropoda melanoleuca*). *PLoS One.* 2015 Sep 16;10(9):e0137712.
23. Yang S, Han Y, Zheng S, Kou X, Asico LD, Huang H, Gao Z, Jose PA, Zeng C. Enhanced Natriuresis and Diuresis in Wistar Rats Caused by the Costimulation of Renal Dopamine D3 and Angiotensin II Type 2 Receptors. *Am J Hypertens.* 2015 Oct;28(10):1267-76.
24. Li F, Feng J, Gao D, Wang J, Song C, Wei S, Qiao M. Shuyu Capsules Relieve Premenstrual Syndrome Depression by Reducing 5-HT3AR and 5-HT3BRExpression in the Rat Brain. *Neural Plast.* 2016;2016:7950781.

Version 2016.10.10